


Pantalla de cobre chapada en estao de 4 ncleos aislada con PFA, cable / cable del sensor PRT



10M



4 CORE
XE-2734-010

25M



4 CORE
XE-2734-025

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pantalla de cobre chapada en estao de 4 ncleos aislada con PFA, cable / cable del sensor PRT

Alta temperatura hasta 260C, construccin de tefln extruido con PFA. Comnmente utilizado para termmetros Pt100, tambin adecuado para otros tipos como Pt130, Pt500 y Pt1000

Especificaciones**Specifications**

General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10m, 25m, 50m, 100m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	4
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket