




Pantalla de cobre chapada en estao de 4 nucleos aislada en PVC aislada, cable / cable sensor PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Pantalla de cobre chapada en estao de 4 ncleos aislada en PVC aislada, cable / cable sensor PRT

Cable aislado de PVC para extender termómetros de resistencia al platino. Comnmente utilizado para termómetros Pt100, pero también adecuado para otros tipos como Pt130, Pt500 y Pt1000

Especificaciones**Specifications**

General Description	4 core PVC insulated 7/0.2mm tin plated copper conductors, twisted cores, polyester tape & nylon cord, tin plated copper braided screen with PVC insulated outer jacket
Cable Length	10, 25, 50, 100m
Cable Type	Mylar shield with tin plated copper braiding
Core / Strands	7/0.2mm
Insulation Material	PVC insulation
Max. Temperature	+105°C
Min. Temperature	-10°C
Number of Wires	4
Standards Met	Cores colour coded 2 x Red / 2 x White to IEC-751 with black outer jacket