

Accesorios de compresin de acero inoxidable - Rosca cnica (BSPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores minerales aislados y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de acero inoxidable con roscas BSPT (British Standard Pipe Taper) se utilizan ampliamente para conectar tuberías y tuberías en diversos entornos exigentes.

Las roscas BSPT son cónicas, lo que significa que el diámetro de la rosca disminuye hacia el extremo, lo que mejora el sellado al apretar el accesorio.

El uso del acero inoxidable proporciona una excelente resistencia a la corrosión, durabilidad y resistencia, lo que hace que estos accesorios sean ideales para condiciones adversas que se encuentran en industrias como el procesamiento químico, el petróleo y gas, y aplicaciones marinas.

Especificaciones

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	Stainless steel, BSPT tapered thread
Thread	1/8"BSPT, 1/4"BSPT, 1/2"BSPT, 3/8"BSPT or 3/4"BSPT available
Probe Diameter	Suitable for both metric and imperial probe diameters from 1.0mm to 12.0mm