

## TermôfÂ³stelos perforados sÃfÂ³lidos



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermopockets de acero inoxidable de varios tamanos están diseñados para su uso con sondas de temperatura de estilo industrial para proporcionar protección para la cubierta de medios corrosivos y para facilitar el reemplazo de la sonda sin perturbar el proceso.

A partir de stock, las conexiones de proceso son BSPP de 1/2 pulgada de serie, se requiere un accesorio de compresión BSPP de 1/8 pulgadas para adaptarse al diámetro de la sonda para conectar el sensor al termobolsillo.

## Specifications

<b>General Description</b>	From stock, process connections are 1/2 inch BSPP as standard, a 1/8 inch BSPP compression fitting to suit diameter of probe is required to connect sensor to thermopocket.
<b>Material</b>	Stainless steel
<b>Dimensions</b>	303SS, 12mm Stem
<b>Type &amp; Stem Diameter</b>	303SS, 12mm Stem
<b>Stem Length</b>	50mm, 100mm, 150mm
<b>Probe Diameter</b>	50mm, 100mm, 150mm stem length