


Termopar de unin expuesta IEC aislada con doble torsione con PTFE y tapn miniatura - Tipo K - 2 M de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar tipo K, J o T. Unin expuesta soldada de respuesta rpida disponible en la especificacin de cable Clase 1 conforme a la BS EN 60584-3:2008 con terminacin de tapn miniatura.

Ideal para aplicaciones de prueba y desarrollo.

Especificaciones

Specifications

| | |
|----------------------------|--|
| Product Code | XE-3506-001 |
| General Description | Fast response, welded exposed junction |
| Sensor Type | Exposed Junction |
| Thermocouple Type | K IEC |
| Cable Length | 2M |
| Cable Type | Teflon® PTFE insulated twin-twisted lead |
| Core / Strands | 1/0.2mm |
| Termination Type | Miniature flat pin plug colour coded in accordance IEC 584 |
| Max. Temperature | 260°C |
| Min. Temperature | -75°C |
| Tolerance | Class 1 |
| Standards Met | IEC-584-3 |
| Connector Type | Plug |
| Connector Size | Miniature |
| Series | Exposed Junctions |