



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar mineral aislado IEC con tapn estndar - Tipo K - 1,5 mm x 150 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Necesitas calibrar tu sensor de temperatura?**

Rastreable a estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de calibración trazable en 3 o 5 puntos utilizando nuestra instalación de calibración interna.

Ten en cuenta que si añades una opción de calibración a tu pedido, te contactaremos para información de entrega.

**Especificaciones**

**Specifications**

<b>Product Code</b>	XE-3700-001
<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral Insulated
<b>Thermocouple Type</b>	K
<b>Element/hot junction</b>	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
<b>Probe Diameter</b>	1.5mm
<b>Probe Length</b>	150mm
<b>Sheath Material</b>	Stainless steel 310
<b>Termination Type</b>	Standard round pin plug, colour coded in accordance IEC 584-3
<b>Max. Temperature</b>	1100°C
<b>Min. Temperature</b>	-40°C
<b>Standards Met</b>	In accordance with IEC 584-1&2
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Colour	Green
--------	-------