

**Termopar mineral aislado de estilo industrial con
cabeza terminal KNE - Tipo K - 6,0 mm de diametro
- 400 mm de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Necesitas calibrar tu sensor de temperatura?

Rastreable segn los estndares de calibracin ISO17025 (UKAS), Labfacility ofrece la opcin de calibracin trazable de 3 o 5 puntos utilizando nuestra instalacin de calibracin interna; puedes seleccionar cualquier punto de temperatura entre -10C y +200C.

Ten en cuenta que si aades una opcin de calibracin a tu pedido, te contactarn para informacin de entrega.

Especificaciones**Specifications**

Product Code	XE-3254-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral insulated Thermocouple fitted with KNE terminal head
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	400mm
Sheath Material	310 stainless steel sheath
Head Type	KNE, Cast aluminium finish, screw type lid with chain
Gland	M16 x 1.5mm Plated brass gland fitted
Max. Temperature	1100°C
Min. Temperature	-40°C
IP Rating	IP67