

**Termopar mineral aislado de estilo industrial con  
cabeza terminal KNE - Tipo K - 6,0 mm de diametro  
- 600 mm de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Necesitas calibrar tu sensor de temperatura?**

Rastreable según los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de calibración trazable de 3 o 5 puntos utilizando nuestra instalación de calibración interna; puedes seleccionar cualquier punto de temperatura entre -10C y +200C.

Ten en cuenta que si añades una opción de calibración a tu pedido, te contactaremos para información de entrega.

**Especificaciones****Specifications**

<b>Product Code</b>	XE-3255-001
<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral insulated Thermocouple fitted with KNE terminal head
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Insulated hot junction
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	600mm
<b>Sheath Material</b>	310 stainless steel sheath
<b>Head Type</b>	KNE, Cast aluminium finish, screw type lid with chain
<b>Gland</b>	M16 x 1.5mm Plated brass gland fitted
<b>Max. Temperature</b>	1100°C
<b>Min. Temperature</b>	-40°C
<b>IP Rating</b>	IP67