


**Termopar de unin expuesta aislado con PFA ANSI  
con enchufe en miniatura - Tipo K - Tamao de  
hebra de 1/0,3 mm - Longitud de 10 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Product Code</b>        | XE-3572-001                             |
| <b>General Description</b> | Fast response, welded exposed junction  |
| <b>Sensor Type</b>         | Exposed Junction                        |
| <b>Thermocouple Type</b>   | K                                       |
| <b>Cable Length</b>        | 10M                                     |
| <b>Cable Type</b>          | Teflon® PFA insulated twin-twisted lead |
| <b>Core / Strands</b>      | 1/0.3mm                                 |
| <b>Termination Type</b>    | Bare Tails                              |
| <b>Max. Temperature</b>    | 250°C                                   |
| <b>Min. Temperature</b>    | -75°C                                   |
| <b>Standards Met</b>       | ANSI MC96.1                             |