

**Termopar de unin expuesta aislado con ANSI PFA
con enchufe en miniatura - Tipo T - Tamao de
hebra de 1 / 0.3 mm - Longitud de 10 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3574-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	T
Cable Length	10M
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI MC96.1
Series	Exposed Junctions