




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de unin expuesta aislado con PFA ANSI con colas desnudas - Tipo K - Tamao de hebra de 1 / 0.3 mm - Longitud de 1 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3543-001
General Description	Fast response, welded exposed junction
Sensor Type	K
Thermocouple Type	Exposed Junction
Cable Length	1M
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.3mm
Termination Type	Bare Tails
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI MC96.1
Series	Exposed Junctions