




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Termopar de unin expuesta aislado con fibra de vidrio ANSI con colas desnudas - Tipo K - 5M de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                              |                                                      |
|------------------------------|------------------------------------------------------|
| <b>Product Code</b>          | XE-3556-001                                          |
| <b>General Description</b>   | Fast response, welded exposed junction               |
| <b>Sensor Type</b>           | Good temperature resistance (not suitable in fluids) |
| <b>Thermocouple Type</b>     | K                                                    |
| <b>Thermocouple Junction</b> | Exposed                                              |
| <b>Cable Length</b>          | 5M                                                   |
| <b>Cable Type</b>            | glassfibre insulated flat pair cable                 |
| <b>Core / Strands</b>        | 1/0.3mm                                              |
| <b>Cable Termination</b>     | Bare tails                                           |
| <b>Max. Temperature</b>      | +350°C (short periods up to 400°C)                   |
| <b>Min. Temperature</b>      | -60°C                                                |
| <b>Tolerance</b>             | Class 2                                              |
| <b>Standards Met</b>         | ANSI MC96.1                                          |