


Detectores de hilo enrollados de 100 Ohmios - Clase A - 2,8 mm x 15 mm



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Por favor, ten en cuenta: los plazos de entrega de los DIN de la clase 1/10 variarn; por favor, consulta o sers contactado con un plazo tras realizar tu pedido.

Especificaciones**Specifications**

Product Code	XE-4746-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class A
Dimensions	2.8mm Width x 15mm Length
Series	RTD Detectors