

Detectores de bobinado de alambre de 100 ohmios - Clase A - 2,8 mm x 25 mm



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Tenga en cuenta: Los plazos de entrega de DIN de clase 1/10 variarn: consulte o se le contactar con un plazo de entrega despues de realizar su pedido.

Caractersticas tcnicas

Specifications

Product Code	XE-4747-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class A
Dimensions	2.8mm Width x 25mm Length
Series	RTD Detectors