

## **Detectores de hilo enrollados de 100 Ohmios - Clase 1/10 DIN - 1,5 mm x 8 mm**



100 OHM WIRE WOUND DETECTOR



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Por favor, tenga en cuenta: los plazos de entrega de DIN de la Clase 1/10 variarn; por favor, consulte para que podamos informarnos antes de realizar su pedido.

**Especificaciones****Specifications**

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Product Code</b>        | XE-4788-001                                                  |
| <b>General Description</b> | Wirewound construction for optimum performance and stability |
| <b>Sensor Type</b>         | Pt100 Single Element                                         |
| <b>Max. Temperature</b>    | 650°C                                                        |
| <b>Min. Temperature</b>    | -200°C                                                       |
| <b>Tolerance</b>           | 1/10DIN relative to IEC 60751 over the range 0°C to 150°C    |
| <b>Dimensions</b>          | 1.5mm Width x 8mm Length                                     |