




Detectores de doble elemento enrollado de 100 Ohmios - Dplex clase B - 1,5 mm x 15 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Especificaciones**Specifications**

Product Code	XE-4810-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Dual element Pt100
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class B to IEC 60751
Dimensions	1.5mm Width x 15mm Length
Series	RTD Detectors