


Detectores de alambre de doble elemento de 100 ohmios - Clase B Dplex - 2.8mm x 25mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-4817-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Dual element Pt100
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class B to IEC 60751
Dimensions	2.8mm Width x 25mm Length
Series	RTD Detectors