




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Detectores de alambre de doble elemento de 100 ohmios - Clase A duplex - 2.8mm x 25mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caractersticas tcnicas****Specifications**

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Product Code</b>        | XE-4847-001                                                  |
| <b>General Description</b> | Wirewound construction for optimum performance and stability |
| <b>Sensor Type</b>         | Dual element Pt100                                           |
| <b>Max. Temperature</b>    | 650°C                                                        |
| <b>Min. Temperature</b>    | -200°C                                                       |
| <b>Tolerance</b>           | Class A to IEC 60751                                         |
| <b>Dimensions</b>          | 2.8mm Width x 25mm Length                                    |
| <b>Series</b>              | RTD Detectors                                                |