

Sonda de uso general con aislamiento mineral de mano - Tipo K,T,PRT



[Click here to view product](#)

Product Code	XE-3806
General Description	Rounded tip suitable for semi solid and liquid temperature measurement
Thermocouple Type	K or T
Calibration	IEC, ANSI or JIS Calibration or PRT (Black)
Probe Diameter	1.5mm, 3.0mm & 6.0mm
Probe Length	100mm, 120mm & 300mm
Cable Length	2m Coiled Cable
Termination Type	Fitted Mini Plug
Max. Temperature	+850°C