


Termopar de Banda Magnética con Enchufe en Miniatura Tipo K IEC - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un termopar de banda magnética flexible capaz de unirse a la mayoría de las superficies ferromagnéticas, la gran resistencia del imán garantiza una buena adhesión, lo que permite una instalación simple. útil en una amplia gama de aplicaciones de monitoreo de temperatura, particularmente en pruebas y desarrollo.

Características técnicas

Specifications

Product Code	XE-3455-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	2M
Cable Type	PTFE insulated wires, twin twisted
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	100°C
Min. Temperature	-50°C
Standards Met	IEC
Series	Magnet Thermocouple