

**Termopar estilo Jubilee Clip de superficie de tubera IEC, trenza de acero inoxidable de fibra de vidrio - Tipo K - para adaptarse a un diametro de tubera de 30-45 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                              |                                                                                 |
|------------------------------|---------------------------------------------------------------------------------|
| <b>Product Code</b>          | XE-3916-001                                                                     |
| <b>General Description</b>   | Stainless steel pipe clamp sensor, jubilee clip style, with adjustable diameter |
| <b>Thermocouple Type</b>     | K                                                                               |
| <b>Thermocouple Junction</b> | Grounded tip for fast response                                                  |
| <b>Cable Length</b>          | 2M                                                                              |
| <b>Cable Type</b>            | Glassfibre stainless steel braided extension cable                              |
| <b>Core / Strands</b>        | 7/0.2mm                                                                         |
| <b>Termination Type</b>      | Bare Tails                                                                      |
| <b>Max. Temperature</b>      | 350°C                                                                           |
| <b>Standards Met</b>         | IEC                                                                             |
| <b>To Fit Pipe Diameters</b> | 30-45mm                                                                         |