

Termopar de punta hipodrmica IEC, PTFE - Tipo T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-9511-001
Sensor Type	chisel point Thermocouple
Thermocouple Type	T
Thermocouple Junction	Grounded hot junction
Probe Diameter	1.6mm
Probe Length	100mm
Sheath Material	316 stainless steel sheath
Cable Length	1.5m
Cable Type	PTFE insulated cable
Core / Strands	7/0.2mm
Termination Type	Bare wire tails (+Brown/-White)
Max. Temperature	+250°C
Min. Temperature	-75°C
Series	Hypodermic Tip Thermocouple