

Interruptor selector de termopar IEC de 6 vas, montaje en banco



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un selector de 6 posiciones en una carcasa moldeada con paneles delanteros y traseros de aluminio. Las entradas traseras se conmutan mediante un selector rotatorio montado delantero a la nica toma de salida.

Cableado internamente con cable compensador de termopar y disponible con tipos de termopar K, J, T, N Y R/S.

ESPECIFICACIONES

Construccin - Unidad porttil / de montaje en banco en carcasa moldeada con paneles delanteros y traseros de aluminio

Tamao - 130 mm de ancho x 60 mm de alto x 145 mm de profundidad

Conectores de entrada 6 conectores de termopar compensados en miniatura

Conectores de salida - Un solo tocacielos de termopar compensado en miniatura

Disponible en los tipos CEIC K, J, T, N & R/S

Specifications