

INTERRUPTOR Selector de termopar IEC 12 Way, montaje en banco



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un selector de 12 vías en una carcasa moldeada con paneles delanteros y traseros de aluminio. Las entradas traseras se conmutan a través de un selector giratorio montado en la parte delantera a la toma de salida única.

Cableado interno con cable compensador de termopar y disponible con termopar tipos K, J, T, N y R/S.

Especificación:

Construcción - Unidad de montaje portátil / banco en carcasa moldeada con paneles delanteros y traseros de aluminio

Tamaño - 130mm de ancho x 60mm de alto x 145mm de profundidad

Conectores de entrada – 12 tomas de termopar compensadas en miniatura

Conectores de salida - Toma de termopar compensada en miniatura única

Disponible en IEC tipo K, J, T, N y R/S

Specifications