



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Certificado de conformidad

ISO9001

CERTIFICATE OF CONFORMITY



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un documento certificado por una autoridad competente de que el bien o servicio suministrado cumple con las especificaciones requeridas.

También se llama certificado de conformidad o certificado de cumplimiento.

Specifications

Product Code	PO-COC
--------------	--------