


Termómetro de resistencia higiénico Pt100, ajuste de estilo RJT de 1.5"



[Click here to view product](#)

Product Code	XE-5525
General Description	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
Sensor Type	Pt100, 3 or 4 wire single element class A
Head Type	KNEF
Connection Type	Pt100, 3 or 4 wire single element class A
Transmitter	Optional SEM-203p transmitter, ranged -50°C to +150°C
Support Dia (D2)	8.0mm
Support Length (L2)	50mm