

Pt100 3 cable clase B Termómetro de resistencia, cabeza KNE con transmisor



[Click here to view product](#)

Product Code	XE-5506
General Description	Pt100 probe with KNE terminal head fitted with a transmitter
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 3 wire
Process Connection	1/2"BSP x 6mm
Calibration	IEC
Probe Diameter	6mm
Probe Length	150mm
Sheath Material	316 stainless steel
Head Type	KNE
Max. Temperature	+450°C (connection head @ 170°C)

Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Transmitter	Yes
Range	Various ranges of transmitter available
Number of Wires	3