

Carcasa montada en pared remota de 4-20mA, entrada Pt100 de 1 metro de plomo RTD / PRT



[Click here to view product](#)

Product Code	XE-5520
General Description	Compatible with Pt100 platinum resistance thermometers
Sensor Type	Surface Mounting for remote installtion
Cable Length	1m
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Thread	Entry Threads M16
Transmitter	Yes
Range	-50°C to +150°C or 0°C to +400°C (either can be re-ranged if desired)
Output	4-20 mA output