




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Detectores de resistencia de platino con cables extendidos, cerámica bobinada - 2 hilos - 300 mm de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Rendimiento ptimo, adecuado para aplicaciones de aire e inmersin (adecuado solo en atmferas secas y no contaminantes donde no estn protegidas. Se requiere una proteccin/revestimiento adecuado para la inmersin en fluidos, sustancias corrosivas, procesos, etc.).

Caracterstics tcnicas

Specifications

Product Code	XE-5595-001
General Description	Platinum Resistance Detectors with extended leads, Wire-Wound Ceramic
Sensor Type	Platinum Resistance, Pt100 Ohm
Sheath Material	Platinum coil wire-wound construction sealed inside a high purity alumina ceramic body
Cable Length	300mm
Cable Type	Glass-fibre insulated nickel
Cable Termination	Bare tails
Max. Temperature	350°C
Min. Temperature	-50°C
Standards Met	European / British Standard Conforms to IEC 60751:2008 / BS EN 60751:2008
Accuracy	Class B
Extended lead wires	2 Wire
Detector Dimensions	2.8mm x 15mm
Fundamental interval	38.5 ohms nominal (100 ohms)
Ice point resistance	100 Ohms
Series	Wire Wound Detector with Extended Leads