

Termmetro infrarrojo IR Temp 2



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Uso industrial: precisión no adecuada para uso médico.

Certificado de calibración rastreable gratuito

Registra las temperaturas máximas/mínimas

Muestra la temperatura ambiente

Termómetro infrarrojo RayTemp 2 - Diseado y fabricado aquí en el Reino Unido, el nuevo termómetro infrarrojo sin contacto RayTemp 2. Simplemente apunta el termómetro al objetivo y pulsa el botón de medición para mostrar la temperatura superficial en el rango de -49,9 a 349,9 °C con una resolución de 0,1 °C o 1 °C y una precisión asegurada de 1 °C en el rango de 0 a 100 °C.

Incluye certificado de calibración rastreable GRATUITO.

Especificaciones

Specifications

Product Code	XE-3941-001
General Description	The RayTemp® 2 incorporates a 5:1 optic ratio (target distance/diameter ratio) and a fixed emissivity of 0.95 making it suitable for a wide range of food and industrial applications.
Type	infrared
Accuracy	±1 °C (0 to 100 °C) otherwise ±2 °C or ±2 % of reading, whichever is greater
Resolution	0.1°C & 1°C
Emissivity	0.95 fixed
Dimensions	25 x 56 x 128mm
Input / Ranges	-49.9 to 349.9°C
Weight	130 grams
Power Consumption	5000 hours continuous use
Power Supply	3 x 1.5 volt AAA



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Important Information

Not suitable for medical use. Please note that certificates supplied do not display an expiry date. They state the date the item was manufactured. The frequency of recalibration is determined by each customer's specific application requirements and intern