

Termmetro industrial Therma 1 (Tipo K)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TENGA EN CUENTA: El plazo de entrega ser de 7 a 10 das laborables al elegir la calibracin

El Therma 1 cuenta con pantallas LCD grandes y fciles de leer con circuito abierto 'Err' y indicacin de batera baja. Cada termmetro funciona con tres pilas AAA de 1,5 voltios que ofrecen una duracin mnima de cinco aos. Las unidades se apagan automticamente tras diez minutos, maximizando la duracin de la batera. Esta funcin puede ser desactivada por el usuario, si es necesario.

No incluye ninguna sonda, por favor elige entre la amplia gama de sondas de Labfacility para adaptarse a este producto

Sensores de temperatura porttiles opcionales para este dispositivo estn disponibles en el men desplegable: 26 sondas disponibles, selecciona tu tipo, longitud y dimetro de la sonda Tipo K.

Especificaciones

Specifications

Product Code	XE-3988-001
General Description	Housed in a robust ABS case that contains 'Biomaster' additive that helps to reduce bacterial growth.
Type	Industrial Digital Thermometer
Thermocouple Type	Type K
Accuracy	$\pm 0.4^{\circ}\text{C} \pm 0.1\%$
Input / Ranges	0.1°C -99.9 to 299.9°C. 1°C 300 to 1372°C
Resolution of data display	0.1°C & 1°C
Dimensions	25 x 56 x 128mm
Power Supply	Battery 3 x 1.5 volt AAA. Battery life 10000 hours/min 5 years
Sensor Type	Sensor type K miniature thermocouple
Weight	130 grams



Accessories

Please note that certificates supplied do not display an expiry date. They state the date the item was manufactured. The frequency of recalibration is determined by each customer's specific application requirements and internal quality procedures.