



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de pernos IEC, acero inoxidable de fibra de vidrio sobre cable trenzado - Tipo K - Paso de rosca M8x1mm - Longitud de perno de 13 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-5563-001
General Description	Bolt sensors are generally screwed into the surface which normally increases the speed of response and sensitivity of the measurement. Suitable for use up to +250°C.
Thermocouple Type	K
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	250°C
Standards Met	IEC
Bolt thread/pitch	M8 x 1mm
Series	Bolt Thermocouple