

Termopar de parche de silicona IEC, PFA con enchufe en miniatura - Tipo K - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-5619-001
General Description	Resistant to many oils & chemicals
Sensor Type	Moulded silicone patch thermocouple, 30 x 15mm patch, self adhesive with aluminium foil backing
Thermocouple Type	K
Cable Length	2M
Cable Type	PFA insulated twin-twist wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	150°C
Min. Temperature	-50°C
Standards Met	IEC