


**Termopar de bayoneta ajustable IEC, cable
trenzado de acero inoxidable de fibra de vidrio -
Tipo J - Terminacin de enchufe estndar con
abrazadera de cable - 3M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas

Specifications

Product Code	XE-5636-001
General Description	Thermocouple with adjustable insertion depth
Sensor Type	6.0mm diameter tip
Thermocouple Type	J
Thermocouple Junction	Grounded tip for fast response
Cable Length	3M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Cable Termination	Standard Plug
Max. Temperature	350°C
Min. Temperature	-60°C
Standards Met	IEC-584-3
Bayonet cap	11.5mm I.D. single slot
Connector Type	Plug
Connector Size	Standard
Thread	M6 x 1mm thread
Spring	200mm
Series	Hermetically Sealed