


Conector de termopar para montaje de PCB IM-K-PCB Tipo K IEC con pines chapados en estao



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Conector de termopar para montaje de PCB IM-K-PCB Tipo K IEC con pines chapados en estao

Diseado especialmente para PCB, pines slidos chapados en estao para montaje directo en placa y provisin para sensor CJC

Especificaciones**Specifications**

Product Code	XE-1109-TIN
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	K IEC
Connector Type	PCB Socket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
Pins	Tin Plated
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS