



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Termopar doble con aislamiento mineral IEC con enchufe dplex tipo K estndar - 3,0 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Necesita calibrar su sensor de temperatura?

Trazable a ISO17025 estándares de calibración (UKAS), Labfacility ofrece una opción de calibración trazable de 3 o 5 puntos utilizando nuestra instalación de calibración interna.

Tenga en cuenta que si agrega una opción de calibración a su pedido, se comunicará con usted con respecto a la información de entrega.

Características técnicas**Specifications**

Product Code	XE-2012-001
General Description	310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes.
Sensor Type	Mineral insulated Type 'K' with Duplex Thermocouple Plug
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	3.0mm
Probe Length	500mm
Termination Type	Standard duplex plug termination (200°C)
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	Conforms to IEC 584 specification
Connector Type	Plug
Connector Size	Standard