

Sonda Pt100 de uso general con cable PFA - Tipo PRT - Dimetro de la sonda de 4,5 mm - Longitud de la sonda de 75 mm - Longitud del cable de 2 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda Pt100 de uso general con cable PFA - Tipo RTD / PRT

Sonda multipropósito Pt100 según IEC751

Clase B, conexión de 4 hilos

Funda rígida de acero inoxidable 316

Plomo apantallado con aislamiento de teflón

Clasificación de temperatura del sensor: de -50 °C a +250 °C

Características técnicas**Specifications**

Product Code	XE-3617-001
General Description	Multipurpose Pt100 Probe to IEC751
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Sheath Material	4.5mm diameter 316 Stainless Steel Sheath
Cable Length	2M
Cable Type	Flexible Teflon® insulated, screened
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC751