




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Detectores de hilo enrollados de 100 ohmios - Clase 1/10 DIN - 2,8 mm x 15 mm



100 OHM WIRE WOUND CLASS 1/10 DIN



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Por favor, tenga en cuenta: los plazos de entrega de DIN de la Clase 1/10 variarn; por favor, consulte para que podamos informarnos antes de realizar su pedido.

Caractersticas tcnicas**Specifications**

Product Code	XE-4796-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	1/10DIN relative to IEC 60751 over the range 0°C to 150°C
Dimensions	2.8mm Width x 15mm Length