

Termopares de horno de cerámica para alta temperatura (1400 oC) con vainas de porcelana aluminosas - Tipos K,N



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

El tubo de soporte de acero inoxidable de 100 mm de largo fija la cubierta de cerámica al cabezal industrial y se incluye en la longitud total de la vaina.

Specifications

General Description	Ceramic thermocouples for furnaces, kilns and other high temperature applications
Sensor Type	Made with rugged 20 gauge (0.8mm dia) thermocouple wire for long-term performance
Thermocouple Type	Types K and N with Aluminous Porcelain sheaths, rated to 1400°C
Sheath Material	12 mm diameter ceramic sheath
Max. Temperature	+1400°C
IP Rating	IP67 rated industrial protection head with cable entry and screw terminal block
Support Tube	21 mm O.D. support tube, 100 mm long