




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

BS (British Standard) Termopar de Disco de Alambre Fino, PTFE con Enchufe en Miniatura - Tipo K - 1M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-5642-001
General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	1M
Cable Type	1/0.2mm PTFE insulated twisted wire
Termination Type	Miniature flat pin plug, colour coded
Standards Met	BS EN 60584-1:2013