

Portátil entre el paquete / sonda de alimentos planos - Tipo K,T



[Click here to view product](#)

Product Code	XE-3848
Sensor Type	Hand Held Between Pack/Flat Food Probe
Thermocouple Type	K or T
Calibration	IEC
Thermocouple Junction	Single element, grounded junction offering fast response.
Probe Diameter	4mm
Probe Length	140mm approximately below the handle
Cable Length	2 metre (extended)
Cable Type	Polyurethane
Cable Termination	moulded on miniature plug.

Max. Temperature	+260°C
Min. Temperature	-50°C
Standards Met	IEC 584