



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Termopar de clip de cocodrilo IEC con cable trenzado de acero inoxidable de fibra de vidrio - Tipo K - 1M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caractersticas tcnicas**Specifications**

Product Code	XE-5652-001
General Description	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	1M
Cable Type	Stainless steel over-braided Glassfibre
Cable Termination	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC
Series	Crocodile Clip Thermocouple