


**Termopar de unin expuesta IEC de par plano
aislado PFA con tapn miniatura - Tipo K - Tamao
de hilo 1/0,3 mm - longitud 3 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Buena resistencia y flexibilidad mecánica, resistente a aceites, ácidos y otros fluidos adversos.

Especificaciones

Specifications

Product Code	XE-3486-001
Sensor Type	Fast response, welded exposed junction
Thermocouple Type	K IEC
Cable Connection	Terminated with a miniature flat pin plug
Cable Length	3M
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Standards Met	(cores & jacket) to IEC-584-3
Insulation Material	Insulation rating: -75°C to 260°C
Series	Exposed Junctions