

Jefe de Terminal de KNE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Cabeza estndar KNE, aluminio pintado en plata con tapn de rosca y tornillo de retencin, adecuado para su uso tanto con bloques terminales cermicos como con transmisores.

Especificaciones**Specifications**

Head Type	KNE Standard Head
Material	Aluminium silver painted
Cable Entry	M20 x 1.5mm
Fixing (Centres)	33mm & 40mm
Max. Temperature	170°C
IP Rating	IP67 rated
Process Entry	1/2" BSPP, 3/4" BSPP or M24
Cap Retaining Type	Screw cap and retaining screw