

**Termopar de unin expuesta de PTFE de un solo
disparo IEC con enchufe en miniatura moldeado -
Tipo K - Longitud de 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un termopar aislado con PTFE de respuesta muy rpida, hecho de una construccin nica de cable de un solo disparo que proporciona un aislamiento hermtico al gas, al vapor y al agua a lo largo del cable.

Este sensor es ideal para el monitoreo de temperatura de piezas de prueba en aplicaciones aeroespaciales y de autoclave.

Caractersticas tcnicas

Specifications

Product Code	XE-0030-001
General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is wildly used in autoclaves.
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	2M
Cable Type	PTFE
Core / Strands	1/0.376mm
Termination Type	Terminated with a miniature flat pin plug (Plug max temp 220°C)
Max. Temperature	250°C
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature