



## Conector de termopar miniatura IM-K-M Tipo K IEC



**Miniature in-line Plug**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Conector de termopar miniatura IM-K-M Tipo K IEC**

Con pines planos adaptados al tipo termopar y ideales para uso general, todos los contactos estn polarizados para asegurar la conexin correcta. Estos conectores aceptan cables termopar de hasta 4 mm de dimetro.

Especificaciones

## Specifications

<b>Product Code</b>	XE-1100-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector