


**Termopar magnético de potencia de tracción
(herradura) tipo K IEC de 9 kg, cable aislado de
PFA con trenza de acero inoxidable que termina
en enchufe estándar - 3.0M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3466-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Sensor Type	Type K thermocouple to IEC 584
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Length	3.0M
Cable Termination	PFA Teflon® insulated with stainless steel over-braid, miniature plug termination
Termination Type	Standard Plug
Max. Temperature	+ 250°C
Min. Temperature	-50°C
Standards Met	IEC