


Termopar magnético de botón de extracción IEC tipo K de 0,7 kg, aislado de PTFE con enchufe en miniatura tipo K - 1 m de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3460-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type
Cable Length	1M
Cable Type	PTFE insulated twin-twist wires
Cable Termination	Miniature Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Series	Magnet Thermocouple