


**Termopar magnético de botón de extracción IEC tipo K de 0,7 kg, aislado de PTFE con enchufe estándar tipo K - 1 m de longitud**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Características técnicas****Specifications**

|                              |                                             |
|------------------------------|---------------------------------------------|
| <b>Product Code</b>          | XE-3461-001                                 |
| <b>General Description</b>   | Ideal for ferrous metal surface measurement |
| <b>Thermocouple Type</b>     | K                                           |
| <b>Thermocouple Junction</b> | Grounded at tip                             |
| <b>Magnet</b>                | 1/2" diameter button type                   |
| <b>Cable Length</b>          | 1M                                          |
| <b>Cable Type</b>            | PTFE insulated twin-twist wires             |
| <b>Cable Termination</b>     | Standard Fitted Plug                        |
| <b>Max. Temperature</b>      | +250°C (plug 200°C)                         |