


**Termopar de imn de traccin IEC tipo K de 3 kg,
cable aislado con PFA con trenza de acero
inoxidable que termina en enchufe estndar - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3463-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	1" diameter (3kg pull)
Cable Length	1M
Cable Type	Stainless steel over-braided PFA Teflon®
Cable Termination	Miniature Plug
Max. Temperature	+250°C (200°C Plug)
Standards Met	IEC
Series	Magnet Thermocouple